
Streamlined PACS (Picture Archiving and Communication System)

## Web-enabled Mini-PACS, QXLink

QXLink allows users to access diagnostic images stored in hospitals anytime, anywhere via the web. Certified by FDA and CE, QXLink is truly a reliable and economical PACS solution ideal for small clinics and hospitals.



[QXLink is practical and convenient PACS for any diagnostic images.]

# Applications

QXLink has 60+ measurement tools for Chiropractic, Orthopedic, Podiatry, and Veterinary. Use full measurement tools such as George's line, cobb angle, spondylolisthesis, cervical/lumbar curve, sagittal spine alignment, and many more! Each tool has intuitive graphic guide for supporting easy diagnosis.



[Chiropractic]

[Orthopedic]

[Podiatric]

[Veterinary]

### Features

#### **Convenient Functionality**

QXLink makes your diagnostic process much easier and more accurate with a variety of convenient functions such as diverse viewing modes and measurement tools with graphic guides.

- Specific Viewing Modes for CT, MRI and MG
- Image Filters & Pseudo Color
- DICOM SR (Structured Report)



[Human Limb Deformity Analysis with Guide]



• Flexible display layout

• Multi-viewers up to 20

mammography view

• 60+ Measurement Tools with graphic guide for Chiropractic, Orthopedic, Podiatry, and Veterinary



[Veterinary TPLO with Guide]



[Report Preview]

#### Maximum Customization

Users can enjoy customized experience of using up to 4 monitors in UHD resolution and customizing frequently used toolbars. They can also adjust display layouts that best suit their usages.

- 4 monitors up to UHD resolution
- Customizable toolbars



[Customizable Toolbars]



[Flexible Display Layout]

#### Remarkable Usability

QXLink offers remarkable usability by supporting hanging protocol, multiple languages, up to 20 viewers, and powerful and diverse view modes.

- Hanging protocol
- Multi-language English, French, Polish, Korean, Chinese

[View Mode: Cross Reference Cine]



• Powerful view modes - series, cine, cross reference,

[View Mode: Mammography]

41–3, Burim–ro 170 beon–gil, Dongan–gu, Anyang–si, Gyeonggi–do, 14055 Republic of Korea Tel +82–70–7011–6161 Fax +82–31–386–8631 E–mail sales@vieworks.com Web xrayimaging.vieworks.com

[Hanging Protocol: One Study List

& Two Viewers]